

## LITTLE PUMP® / LITTLE PUMP® DUAL ELECTROMECHANICAL PNEUMATIC TOURNIQUETS

**LITTLE PUMP®** range, DDM electromechanical tourniquets made in France !



**LITTLE PUMP®**



**LITTLE PUMP® DUAL**

- Ease of use
- French quality

CE 0459

## LITTLE PUMP® TOURNIQUET

DDM **LITTLE PUMP®** is an electromechanical pneumatic tourniquet providing a bloodless surgical field for upper and lower limb surgery in complete safety.

**LITTLE PUMP®** is available in two models:

- **LITTLE PUMP®** : with 1 pressure circuit for using with one cuff
- **LITTLE PUMP® DUAL**: with 2 independent pressure circuits for using with two cuffs



### Ease of use

- A start button on the side of the device
- Simple pressure adjustment with + or - keys, then validate
- Simple timer setting
- Digital pressure and timer display
- Color coding to differentiate the two circuits (for LITTLE PUMP® DUAL)
- Programmable audible and visual alarms
- ALRIV function



### Made in France

All our tourniquets are entirely developed and manufactured in France at our production site in Agen.

## REFERENCES

REFERENCE	DESCRIPTION
G10705	Little Pump® electro-mechanical tourniquet with 1 regulated pressure circuit
G10706	Little Pump® Dual electro-mechanical tourniquet with 2 regulated pressure circuits

## AN ADAPTED CUFF RANGE FOR ALL TYPES OF PATIENTS

### Monobloc PU arm and lower limb cuffs

Monobloc PU arm and lower limb cuffs, decontaminable with 1 or 2 bladders, made of a single Monobloc unit with a ribbon to secure the Velcro or a fastening strap (optional) and featuring a range of sizes identifiable by color-coding of the tubing.

Available in straight and conical design to better fit the limb.



	WITH STRAP		WITH STRAP		WITH STRAP	
	1 BLADDER - STRAIGHT		1 BLADDER - CONICAL		2 BLADDERS - STRAIGHT	
Infant arm cuff	GBDM101	GBDSM101	/	/	/	/
Child arm cuff	GBDM102	GBDMS102	/	/	GBDM202	GBDMS202
Small adult arm cuff	GBDM103	GBDMS103	/	/	/	/
Adult arm cuff	GBDM104	GBDMS104			GBDM204	GBDMS204
Large adult arm cuff	/	/	GBCM105	GBCMS105	/	/
					2 BLADDERS - CONICAL	
Small adult lower limb cuff	GCDM103	GCDMS103	/	/	/	/
Adult lower limb cuff	GCDM104	GCDMS104	GCCM104	GCCMS104	GCCM204	GCCMS204
Large adult lower limb cuff	GCDM105	GCDMS105	GCCM105	GCCMS105	/	/
XL adult lower limb cuff	GCDM106	GCDMS106	GCCM106	GCCMS106	/	/
XXL adult lower limb cuff	GCDM107	GCDMS107	GCCM107	GCCMS107	/	/

### Autoclavable arm and lower limb polyester cuffs

Disinfectable and autoclavable polyester arm and lower limb cuffs are made of a removable silicone pouch for easy autoclaving.

Available in straight and conical design to better fit the limb.



	1 BLADDER - STRAIGHT	2 BLADDERS - STRAIGHT	LOWER LIMB CUFFS - CONICAL
Infant arm/lower limb cuff	GBS101	/	/
Child arm cuff	GBS102	GBS202	/
Small adult arm cuff	GBS103	/	/
Adult arm cuff	GBS104	GBS204	/
Child lower limb cuff	GCS102	/	/
Small adult lower limb cuff	GCS103	/	/
Adult lower limb cuff	GCS104	GCS204	GCCS104
Large adult lower limb cuff	GCS105	/	GCCS105
XL adult lower limb cuff	GCS106	/	GCCS106

## ACCESSORIES










### Wheeled base and basket *ref A10701*

- An arm equipped with a handle for easy movement.
- One foot with wheels and brake.
- A storage basket.



### Connectors

#### CLIP CONNECTION SYSTEM WITH O-RING

REFERENCE	IMAGE	MATERIAL	DESIGNATION
A11012M		METAL	4.8 mm fluted in-line female coupler
A11016		ACETAL	4.8 mm fluted in-line female end-piece with valve
A11008M		METAL EPDM	4.8mm metal male connector with O-ring
A11010		ACETAL EPDM	4.8 mm fluted in-line male end-piece with sealing valve
A11322		EPDM	O-ring for male connector
A20001		METAL	3 m blue extension cable + male and female connectors
A20006		METAL	3 m red extension cable + male and female connectors
A20003		SILICONE + ACETAL	3 m blue autoclavable, underwater extension cable + male and female connectors
A20005		SILICONE + ACETAL	3 m red autoclavable, underwater extension cable + male and female connectors

### DDM FRENCH QUALITY

#### CE mark 0459

Tourniquets are Class IIa medical devices. They comply with current European regulations.

#### After-Sales Service

GSH provides tourniquet maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

#### Warranty

Tourniquets are guaranteed for 2 years (under normal conditions of use). The warranty does not cover faulty maintenance by the user.